



## Level Detection at Drip Chambers

In order to guarantee the correct dosage of liquid quantities during the administration of **infusions** at any time, it is necessary to monitor permanently the fluid level in drip chambers. Based on our years of experience with the special requirements of our medical device customers, we have developed the **drip chamber sensor SONOCHECK ALD**. The detector monitors the liquid level in drip chambers and is, for example, used in dialysis machines or parenteral nutrition. As a **clamp-on** component, the drip chamber is simply inserted into the sensor which can be installed freely suspended or integrated permanently in a device.

The **non-invasive ultrasonic sensor** is characterized by the following features:

Meets the strict safety standards for medical electrical devices according to **IEC 60601-1**

Dimensions and design adaptable to every custom-specific drip chamber

### Detailed Product Information



#### **SONOCHECK ALD**

clamp-on sensor for level detection at drip chambers

---

WORLDWIDE

---

Do you have any questions?

Melissa Öser

+49 (0)345 / 133 17-840

m.oeser@sonotec.de

## THE AMERICAS

---

Do you have any questions?

+1 631 / 415 4758

sales@sonotecusa.com

## NEWS:

---

18.12.2018

**SONOTEC end-of-the-year review: strengthen international networks - foster innovations - invest in the future**

[more News...](#)

## EVENT CALENDAR:

---

05.02.2019

**MD&M West**

[all Events...](#)

## CUSTOMER REVIEWS:

---

**Our customers on our sensors for fluid monitoring**

[all customer reviews...](#)